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CIMMO® HD (fast setting, gray powder)

INTRODUCTION

CIMMO® HD is a dental cement composed of mineral oxides in the form of fine hydrophilic particles. It is indicated in cases of root perforation (endodontic and furcation), root perforation due to internal resorption, retrofilling, direct pulp protection, filling of cavities with enamel without dentin support (dentin substitute).

CIMMO® HD offers the following benefits:

- Particle size that allows complete hydration during mixing;
- Excellent marginal sealing ability that prevents the migration of fluids into the root canal;
- Excellent sealing of root perforations (canal and furcation) by inducing the formation of periradicular cementum.
- Induction of dentin barrier formation when applied over pulp exposures.
- Unlike other cements that require a completely dry field, **CIMMO® HD** is indicated even in places without adequate humidity control, such as in surgeries for (treatment of perforations or retrofills), without losing its properties.

COMPOSITION:

Calcium Oxide, Calcium Carbonate, Magnesium Oxide, Dicalcium Silicate, Aluminum Oxide, Sodium Oxide, Potassium Oxide and Pozzolana (giving rise to calcite granules responsible for remineralizing the region).

PROPERTIES:

- 1 - Initial setting time of approximately 5 minutes and final setting time of 10 minutes;
- 2 - Alkalinity: After spatulation with liquid, it has a pH value of 10, which in 3 hours stabilizes at a value of 12;
- 3 - Radiopacity: **CIMMO® HD** presents isodensity (radiopacity) relative to that of dentin;
- 4 - Resistance to compression: 40 MPa after 24 hours and 65 MPa after 21 days.

TECHNIQUE USED: Sterilize the glass plate, spatula and instruments for insertion and condensation of **CIMMO® HP Cement**.

FOLLOW INSTRUCTIONS FOR EACH USE:

Furcation perforation occurred during canal location:

- 1 - Carry out the suction of the bleeding through an endodontic canal;
- 2 - Divide a measure of the powder in half and place a few drops of the liquid separately. With a small flexible CIMMO spatula, mix to a denser consistency, adding more powder if necessary (Putty consistency). If necessary, wet the Mp of the spatula in the separated drop and moisten the cement so that it becomes hydrated and shiny;
- 3 - Insert the **CIMMO® HD** cement into the perforation with an insertion spatula;
- 4 - Homogenize the surface of the **CIMMO® HD** cement with dry cotton coupled with tweezers. Apply 30 seconds of air from the syringe to speed up the adjustment;
- 5 - The working time should be 3 to 5 minutes;
- 6 - Instrument the canals and perform the endodontic techniques of choice.

Teeth submitted to previous treatments, with previously diagnosed furcation perforation:

- 1 - Cauterize the granulation tissue of the furcation;
- 2 - Divide a measure of the powder in half and place a few drops of the liquid separately. With a small flexible CIMMO spatula, mix to a denser consistency, adding more powder if necessary (Putty consistency). If necessary, wet the Mp of the spatula in the separated drop and moisten the cement so that it becomes hydrated and shiny;
- 3 - Insert the **CIMMO® HD** cement into the perforation with an insertion spatula;
- 4 - Homogenize the surface of the **CIMMO® HD** cement with dry cotton coupled with tweezers. Apply 30 seconds of air from the syringe to speed up the adjustment;
- 5 - The working time should be 3 to 5 minutes;
- 6 - Instrument the canals and perform the endodontic techniques of choice.

Cervical perforations:

- 1- Perform flap detachment;
- 2- In a glass plate, divide a measure of the powder in half and insert a few drops of the separate liquid. With a small flexible CIMMO spatula, combine half of the portion with a drop until reaching a denser consistency, adding more powder if necessary (Putty consistency). If necessary, wet the tip of the spatula in the separated drop and moisten the cement so that it stays hydrated bright;
- 3- Insert the CIMMO® HD cement into the perforation with an insertion spatula;
- 4- Homogenize the surface of the CIMMO® HD cement with dry cotton coupled with clinical tweezers. Apply 30 seconds of air from the triple syringe to accelerate setting;
- 5- The working time should be 3 to 5 minutes;
- 6- Repositioning of flap and suture after 10 minutes of setting.

Retrooturation in Apical Surgery:

- 1- After retro-instrumentation;
- 2- In a glass plate, divide a measure of the powder in half and insert a few drops of the separate liquid. With a small flexible CIMMO spatula, combine half of the portion with a drop until reaching a denser consistency, adding more powder if necessary (Putty consistency). If necessary, wet the tip of the spatula in the separated drop and moisten the cement so that it stays hydrated bright;
- 3- Insert the CIMMO® HD cement in the apical cavity with suitable instruments;
- 4- Homogenize the surface of the CIMMO® HD cement with dry cotton coupled with clinical tweezers. Apply 30 seconds of air from the triple syringe to accelerate setting;
- 5- The working time should be 3 to 5 minutes;

WARNINGS | CONTRAINDICATIONS

Wear goggles, mask and gloves when handling CIMMO® HD Cement. In case of contact with eyes or skin, wash with water. Only open the blister or vial immediately before use.

ATTENTION: The information contained in this package insert is based on scientific clinical and laboratory studies. However, the success of procedures with CIMMO® HD cement depends on a correct diagnosis, careful surgical technique, conditions of the tooth being treated and the patient's systemic condition. This product must be used in accordance with the instructions in this package insert. The CIMMO® HD cement has no adhesiveness in intra-root cores, therefore their cementation must be done with conventional cement indicated for each type of core.

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