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CIMMO® DTA (quick setting, off white powder)

INTRODUCTION

CIMMO® DTA is a dental reparative cement composed of mineral oxides in the form of fine hydrophilic particles. It is indicated as a base for deep restorations, dentin desensitizer, filling of cavities with unsupported enamel (substitute for dentin strengthening the pulp chamber), root framework for installation of fiberglass or metal core.

CIMMO® DTA offers the following benefits:

- Particle size that allows complete hydration during mixing;
- Excellent marginal sealing ability that prevents fluid migration;
- Excellent suitability of the medium allowing rapid resolution in teeth destroyed by caries;
- Unlike other cements that require a completely dry field, CIMMO® DTA cement is indicated in places without adequate control, without losing its properties.

COMPOSITION:

Calcium Oxide, Calcium Carbonate, Magnesium Oxide, Dicalcium Silicate, Aluminum Oxide, Sodium Oxide, Potassium Oxide and Pozzolana (giving rise to calcite granules responsible for remineralizing the region).

PROPERTIES:

- 1- Initial setting time of approximately 5 minutes and final setting time of 10 minutes;
- 2- Alkalinity: After spatulation with liquid, it has a pH value of 10 which in 3 hours stabilizes at a value of 12;
- 3- Radiopacity: CIMMO® DTA cement has low radiopacity and its compressive strength: 40 MPa after 24 hours and 65 MPa after 21 days;
- 4- When mixing the CIMMO® DTA cement, a gray hue can be observed that decreases with the beginning of the setting time, stabilizing in an ice color without risk of darkening of the teeth.

TECHNIQUE USED: Sterilize the glass plate, the spatula and the instruments for inserting and condensing the CIMMO® DTA cement.

FOLLOW INSTRUCTIONS FOR EACH USE:

Base for deep restorations; dentin desensitizer:

- 1- Divide a measure of the powder in half and place a few drops of the liquid separately. With a small flexible CIMMO spatula, mix to a denser consistency, adding more powder if necessary (Putty consistency);
- 2- Insert the CIMMO® DTA cement with an insertion spatula;
- 3- Homogenize the surface of the CIMMO DTA® cement with dry cotton coupled with clinical tweezers;
- 4- Apply on the CIMMO® DTA cement, air from the triple syringe for 30 seconds;
- 5- The working time should be 3 to 5 minutes;
- 6- Composite resin restoration can be performed after 5 minutes.

Cavity filling with unsupported enamel (dentin substitute strengthening the dental crown):

- 1- Divide a measure of the powder in half and place a few drops of the liquid separately. With a small flexible CIMMO® DTA spatula, mix to a denser consistency, adding more powder if necessary (PuRy consistency);
- 2- Insert the CIMMO® DTA cement with an insertion spatula;
- 3- Homogenize the cement surface with dry cotton coupled to clinical tweezers;
- 4- Apply on the CIMMO® DTA cement, air from the triple syringe for 30 seconds;
- 5- The working time should be between 3 and 5 minutes;
- 6- It can be installed as a single element (only CIMMO® DTA cement) or immediate restoration after 5 minutes;
- 7- When used as a single element, in the second session, it is recommended to lower the CIMMO® DTA cement with a spherical diamond bur and perform a composite resin restoration.

Root framework for installation of fiberglass post or metallic post:

- 1- Divide a measure of the powder in half and place a few drops of the liquid separately. With a small flexible CIMMO® spatula, mix to a denser consistency, adding more powder if necessary (PuRy consistency);
- 2- Insert the CIMMO® DTA cement in the root canal, previously prepared for the placement of the post, using Mespaden or Lentulo drill. Immediately place the post (which can be made of fiberglass) with the sole purpose of modeling the cement inserted in the canal to form a framework, before the cement sets;
- 3- Apply on the CIMMO® DTA cement, air from the triple syringe for 30 seconds;
- 4- After 5 minutes, remove the pin;
- 5- Cement the post with resin cement for glass fiber or zinc phosphate cement for metallic post;
- 6- Finish core filling with composite resin;
- 7- Perform prophetic rehabilitation.

WARNINGS AND CONTRAINDICATIONS:

Wear goggles, mask and gloves when handling CIMMO® DTA Cement. In case of contact with eyes or skin, wash with water. Only open the bottle immediately before use.

ATTENTION: The information contained in this package insert is based on scientific, clinical and laboratory studies. However, the success of procedures with CIMMO® DTA cement depends on a correct diagnosis, careful surgical technique, conditions of the tooth being treated and the patient's systemic condition. This product must be used in accordance with the instructions in this package insert. The CIMMO® DTA cement does not adhere to intra-root cores, so the cementation of the posts must be done, after the CIMMO® DTA cement has set, with conventional cement indicated for each type of post.

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